

Anti-Cingulin Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P2M7

Item Number: CM0021864

Introduction

Rabbit polyclonal antibody targeting human cingulin (CGN). Validated for WB, IHC-P, ELISA. Ideal for tight junction biology and epithelial barrier research. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CGN Rabbit pAb
Remarks/Alias	CGN; cingulin
Species	Human
Gene ID (Human)	57530
Gene ID	57530
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 904-1203 of human CGN (NP_065821.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9P2M7
Protein Weight	137kDa
Shipping	Ice bag