

Anti-Cinp (Cyclin-Dependent Kinase 2-Interacting Protein) Polyclonal Antibody For Wb, Elisa - Q9Bw66

Item Number: CM0020405



Introduction

Rabbit polyclonal antibody against human CINP, validated for WB and ELISA. Detects human and rat CINP. Suitable for DNA replication, cell cycle, and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	CINP Rabbit pAb
Species	Human
GeneID (Human)	51550
GeneID	51550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human CINP (NP_116019.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BW66
Protein Weight	24kDa
Shipping	Ice bag