


Anti-Cirbp Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q14011

Item Number: CM0003635



Introduction

CIRBP rabbit monoclonal antibody validated for WB, IP, IF-P, IHC-P, and ELISA. Reacts with human CIRBP, a cold-inducible RNA-binding protein involved in stress response and cell survival.

[Learn More](#)

Parameter	Value
Product Name	CIRBP Rabbit mAb
Remarks/Alias	CIRP; CIRBP
Species	Human
GeneID (human)	1153
GeneID	1153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CIRBP (Q14011)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q14011
Protein Weight	19kDa
Shipping	Ice bag