


Anti-Cirbp/Cirp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14011

Item Number: CM0028951



Introduction

Rabbit polyclonal antibody against human CIRBP (CIRP), validated for WB, IHC-P, ELISA. Cross-reacts with mouse, rat. Target molecular weight 19 kDa; involved in cold-inducible stress response and RNA binding. For research use in genotoxic stress studies.

[Learn More](#)

Parameter	Description
Product Name	CIRBP Rabbit pAb
Remarks/Alias	CIRP; CIRBP
Species	Human
Gene ID (Human)	1153
Gene ID	1153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CIRBP (NP_001271.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14011
Protein Weight	19kDa
Shipping	Ice bag