


Anti-Cirh1A Rabbit Polyclonal Antibody For Wb, Elisa - Q969X6

Item Number: CM0030082



Introduction

Anti-CIRH1A rabbit polyclonal antibody validated for Western blot (WB) and ELISA, with cross-reactivity to human and mouse. Targets UTP4, a ribosome biogenesis factor involved in 18S rRNA processing. Ideal for nucleolar research and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	CIRH1A Rabbit pAb
Remarks/Alias	NAIC; CIRH1A; CIRHIN; TEX292
Species	Human
GeneID (Human)	84916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 497-686 of human CIRH1A (NP_116219.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q969X6
Protein Weight	77kDa
Shipping	Ice bag