


Anti-Cisd1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nz45

Item Number: CM0018434



Introduction

Rabbit polyclonal antibody against human Cisd1 (mitoNEET), suitable for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying mitochondrial iron-sulfur cluster transfer, redox sensing, and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	CISD1 Rabbit pAb
Remarks/Alias	ZCD1; MDS029; C10orf70; mitoNEET; CISD1
Species	Human
GeneID (Human)	55847
GeneID	55847
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-108 of human CISD1 (NP_060934.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NZ45
Protein Weight	12kDa
Shipping	Ice bag