


Anti-Cisd2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N5K1

Item Number: CM0027393



Introduction

CISD2 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA, cross-reacts with Human, Mouse, Rat. Targets CDGSH iron-sulfur domain-containing protein 2, a key autophagy regulator and lifespan-related protein.

[Learn More](#)

Parameter	Value
Product Name	CISD2 Rabbit pAb
Remarks/Alias	ERIS; WFS2; ZCD2; NAF-1; Miner1; CISD2
Species	Human
GeneID (Human)	493856
GeneID	493856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-135 of human CISD2 (NP_001008389.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5K1
Protein Weight	15kDa
Shipping	Ice bag