


Anti-Cisd2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N5K1

Item Number: CM0017416



Introduction

High-quality rabbit polyclonal antibody targeting human CISD2 (ERIS/Miner1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for autophagy and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	CISD2 Rabbit pAb
Remarks/Alias	ERIS; WFS2; ZCD2; NAF-1; Miner1; CISD2
Species	Human
Gene ID (Human)	493856
Gene ID	493856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-135 of human CISD2 (NP_001008389.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5K1
Protein Weight	15kDa
Shipping	Ice bag