

Anti-Cish Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nse2

Item Number: CM0017761

Introduction

CISH Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat CISH (CIS-1, SOCS) involved in JAK-STAT5 signaling. Sourced from rabbit, recombinant protein immunogen. Ideal for research on cytokine regulation and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	CISH Rabbit pAb
Remarks/Alias	CIS; G18; SOCS; CIS-1; BACTS2; CISH
Species	Human
GeneID (human)	1154
GeneID	1154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-130 of human CISH (NP_659508.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NSE2
Protein Weight	29kDa
Shipping	Ice bag