

Anti-Cish Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Nse2

Item Number: CM0010582

Introduction

CISH rabbit monoclonal antibody validated for WB, IP, ELISA in human samples. Detects the ~29 kDa cytokine-inducible SH2-containing protein, a key regulator of JAK-STAT5 signaling. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	CISH Rabbit mAb
Remarks/Alias	CIS; G18; SOCS; CIS-1; BACTS2
Species	Human
GenelD (human)	1154
GenelD	1154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human CISH (NP_659508.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9NSE2
Protein Weight	29kDa
Shipping	Ice bag