

Anti-Cish Rabbit Polyclonal Antibody For Wb & Elisa - Q9Nse2

Item Number: CM0027348

Introduction

Rabbit polyclonal antibody targeting human CISH, validated for WB and ELISA. Cross-reacts with mouse and rat. CISH regulates cytokine signaling via JAK-STAT pathway inhibition, ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CISH Rabbit pAb
Remarks/Alias	CIS; G18; SOCS; CIS-1; BACT52; CISH
Species	Human
GenelD (Human)	1154
GenelD	1154
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NSE2
Protein Weight	29kDa
Shipping	Ice bag