

Anti-Citrate Synthase Monoclonal Antibody For Wb - 075390

Item Number: CM0012012



Introduction

Rabbit monoclonal antibody against human Citrate Synthase (CS), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Key reagent for mitochondrial metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	Citrate synthetase Rabbit mAb
Remarks/Alias	CS; Citrate synthetase
Species	Human
Gene ID	1431
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 367-466 of human Citrate synthetase (075390).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075390
Protein Weight	52kDa
Shipping	Ice bag