

Anti-Citrate Synthetase Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75390

Item Number: CM0007155



Introduction

Rabbit monoclonal antibody against human citrate synthase. Reacts with human, mouse, and rat. Validated for WB, IHC-P, IF/ICC, and ELISA applications. Ideal for mitochondrial metabolism and TCA cycle research.

[Learn More](#)

Parameter	Value
Product Name	Citrate synthetase Rabbit mAb
Remarks/Alias	CS; citrate synthase; Citrate synthetase
Species	Human
Gene ID (Human)	1431
Gene ID	1431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-145 of human Citrate synthetase (NP_004068.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75390
Protein Weight	52kDa
Shipping	Ice bag