

Anti-Citrate Synthetase Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - 075390

Item Number: CM0007118



Introduction

Recombinant rabbit monoclonal antibody against human Citrate synthetase (CS, 075390). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for Krebs cycle and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Citrate synthetase Rabbit mAb
Remarks/Alias	CS; citrate synthase; Citrate synthetase
Species	Human
GeneID (Human)	1595
GeneID	1431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-145 of human Citrate synthetase (NP_004068.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075390
Protein Weight	52kDa
Shipping	Ice bag