


Anti-Ck14 Rabbit Polyclonal Antibody For Wb, If, Elisa - P02533

Item Number: CM0023593



Introduction

CK14 Rabbit Polyclonal Antibody is validated for WB, IF/ICC, IF-P, ELISA, with reactivity for human, mouse, rat CK14. Suitable for epithelial biology and dermatology research applications.

[Learn More](#)

Parameter	Value
Product Name	CK14 Rabbit pAb
Remarks/Alias	KRT14; CK14; EBS3; EBS4; K14; NFJ; keratin 14
Species	Human
GeneID	3861
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 372-472 of human CK14 (NP_000517.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02533
Protein Weight	51kDa
Shipping	Ice bag