

Anti-Ckap4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q07065

Item Number: CM0025729

Introduction

CKAP4 Rabbit Polyclonal Antibody (pAb) validated for WB, IF/ICC, ELISA applications. Recognizes human, mouse, rat CKAP4, a 66 kDa endoplasmic reticulum membrane protein (SWISS Q07065). Ideal for IVD research.

[Learn More](#)

Parameter	Specification
Product Name	CKAP4 Rabbit pAb
Remarks/Alias	p63; CLIMP63; CLIMP-63; ERGIC-63; CKAP4
Species	Human
Gene ID (Human)	10970
Gene ID	10970
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-400 of human CKAP4 (NP_006816.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07065
Protein Weight	66kDa
Shipping	Ice bag