

Anti-Cklf Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ubr5

Item Number: CM0024584

Introduction

CKLF rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat CKLF. Target is a chemokine-like factor involved in inflammation, chemotaxis, and skeletal muscle regeneration. Ideal for immunology and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	CKLF Rabbit pAb
Remarks/Alias	C32; CKLF1; CKLF2; CKLF3; CKLF4; UCK-1; HSPC224; CKLF
Species	Human
GeneID	51192
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-67 of human CKLF (NP_057410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBR5
Protein Weight	17kDa
Shipping	Ice bag