

Anti-Ckm/Ckb Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P06732

Item Number: CM0010618



Introduction

Rabbit monoclonal antibody against human CKM/CKB, validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for energy metabolism studies.

[Learn More](#)

Parameter	Specification
Product Name	CKM/CKB Rabbit mAb
Remarks/Alias	CKMM; M-CK; CPK-M
Species	Human
Gene ID (Human)	1158
Gene ID	1158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3-381 of human CKM/CKB (NP_001815.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06732
Protein Weight	43 kDa
Shipping	Ice bag