


Anti-Ckmt1A Rabbit Polyclonal Antibody For Wb, Elisa - P12532

Item Number: CM0027190



Introduction

CKMT1A Rabbit Polyclonal Antibody validated for WB and ELISA, reacts with human, mouse, rat. Targets mitochondrial creatine kinase U-type, a key enzyme in energy metabolism, suitable for studies on mitochondrial function and cellular energy homeostasis.

[Learn More](#)

Parameter	Value
Product Name	CKMT1A Rabbit pAb
Remarks/Alias	CKMT1; U-MtCK; mia-CK; CKMT1A
Species	Human
GeneID (Human)	548596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 258-417 of human CKMT1A (NP_001015001.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12532
Protein Weight	47kDa
Shipping	Ice bag