


## Anti-Ckmt1B Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P12532

Item Number: CM0026283



### Introduction

Rabbit polyclonal antibody targeting human CKMT1B (mitochondrial creatine kinase U-type). Validated for WB, IF/ICC, IP, ELISA in human, mouse, rat samples. 47kDa. Ideal for energy metabolism studies.

[Learn More](#)

| Parameter        | Specification                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CKMT1B Rabbit pAb                                                                                                |
| Remarks/Alias    | CKMT; CKMT1; UMTCK; CKMT1B                                                                                       |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 1159                                                                                                             |
| GeneID           | 1159                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human CKMT1B (NP_066270.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, IP, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P12532                                                                                                           |
| Protein Weight   | 47kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |