


Anti-Ckmt2 Rabbit Polyclonal Antibody For Wb And Elisa - P17540

Item Number: CM0023995



Introduction

Anti-CKMT2 rabbit polyclonal antibody validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Targets mitochondrial creatine kinase S-type, essential for energy metabolism in cardiac and skeletal muscle research.

[Learn More](#)

Parameter	Value
Product Name	CKMT2 Rabbit pAb
Remarks/Alias	SMTCK; CKMT2
Species	Human
GeneID (Human)	1160
GeneID	1160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-230 of human CKMT2 (NP_001816.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17540
Protein Weight	48kDa
Shipping	Ice bag