


Anti-Ckmt2 (Smtck) Polyclonal Antibody For Wb And Elisa - P17540

Item Number: CM0025542



Introduction

Rabbit polyclonal antibody targeting human CKMT2 (SMTCK), validated for WB and ELISA. Detects human, mouse, and rat. Use for energy metabolism and mitochondrial research in cardiac and skeletal muscle.

[Learn More](#)

Parameter	Value
Product Name	CKMT2 Rabbit pAb
Remarks/Alias	SMTCK; CKMT2
Species	Human
GeneID (human)	1160
GeneID	1160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-230 of human CKMT2 (NP_001816.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17540
Protein Weight	48kDa
Shipping	Ice bag