


Anti-Cks2 Polyclonal Antibody For Wb, Elisa - P33552

Item Number: CM0029182



Introduction

Rabbit polyclonal antibody against human CKS2 (Cyclin-dependent kinases regulatory subunit 2). Validated for WB and ELISA, cross-reacts with rat. Ideal for cell cycle and division research.

[Learn More](#)

Parameter	Value
Product Name	CKS2 Rabbit pAb
Remarks/Alias	CKSHS2; CKS2
Species	Human
GeneID (Human)	1164
GeneID	1164
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human CKS2 (NP_001818.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P33552
Protein Weight	10kDa
Shipping	Ice bag