


Anti-Clasp1 Polyclonal Antibody For Wb, Elisa - Q7Z460

Item Number: CM0026031



Introduction

Rabbit polyclonal antibody against human CLASP1. Suitable for WB and ELISA; cross-reacts with human, mouse, rat. Target is a microtubule plus-end tracking protein involved in cell migration and mitosis. Ideal for cell biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CLASP1 Rabbit pAb
Remarks/Alias	MAST1; CLASP1
Species	Human
GeneID (Human)	23332
GeneID	23332
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1389-1538 of human CLASP1 (NP_056097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z460
Protein Weight	169kDa
Shipping	Ice bag