

Anti-Clathrin Heavy Chain Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q00610

Item Number: CM0009909



Introduction

Rabbit monoclonal antibody targeting human Clathrin heavy chain (CLTC, Q00610). Validated for WB, IP, IF/ICC, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for coated pit and vesicle trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	Clathrin heavy chain Rabbit mAb
Remarks / Alias	Hc; CHC; CHC17; MRD56; CLH-17; CLTCL2
Species	Human
Gene ID (Human)	1213
Gene ID	1213
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1571-1670 of human Clathrin heavy chain (NP_004850.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00610
Protein Weight	187kDa
Shipping	Ice bag