

Anti-Clathrin Heavy Chain Polyclonal Antibody For Wb And Elisa - Q00610

Item Number: CM0029449



Introduction

Clathrin heavy chain Rabbit Polyclonal Antibody validated for WB, ELISA. Recognizes human, mouse, rat. Immunogen: synthetic peptide from aa 1550-1650. Ideal for coated pit, endocytosis, autophagy, and mitosis research.

[Learn More](#)

Parameter	Details
Product Name	Clathrin heavy chain Rabbit pAb
Remarks/Alias	Hc; CHC; CHC17; MRD56; CLH-17; CLTCL2; Clathrin heavy chain
Species	Human
Gene ID (Human)	1213
Gene ID	1213
Immunogen	Synthetic peptide corresponding to amino acids 1550-1650 of human Clathrin heavy chain (NP_004850.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00610
Protein Weight	192kDa
Shipping	Ice bag