

Anti-Claudin 1 Polyclonal Antibody For Wb & Elisa - 095832

Item Number: CM0020466

Introduction

High-quality rabbit polyclonal antibody against human Claudin 1, validated for WB and ELISA. Ideal for studying tight junction biology, barrier function, and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Claudin 1 Rabbit pAb
Remarks/Alias	CLD1; SEMP1; ILVASC; Claudin 1
Species	Human
GeneID (Human)	9076
GeneID	9076
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 137-211 of human Claudin 1 (NP_066924.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	095832
Protein Weight	23 kDa
Shipping	Ice bag