

Anti-Claudin 1 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - O95832

Item Number: CM0007382

Introduction

KD validated anti-Claudin 1 rabbit monoclonal antibody for WB, IP, IF-P, IHC-P, ELISA. Detects human, mouse, rat CLDN1. Tight junction protein, essential for epithelial barrier and viral entry. Ideal for epithelial biology and virology research.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] Claudin 1 Rabbit mAb
Remarks/Alias	CLD1; SEMP1; ILVASC
Species	Human
Gene ID (Human)	9076
Gene ID	9076
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95832
Protein Weight	23kDa
Shipping	Ice bag