

Anti-Claudin-3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15551

Item Number: CM0006785

Introduction

Claudin-3 Rabbit Monoclonal Antibody detects human Claudin-3 (UniProt O15551) with applications in WB, IHC-P, IF/ICC, and ELISA. Essential for tight junction and cell junction research.

[Learn More](#)

Parameter	Value
Product Name	Claudin-3 Rabbit mAb
Remarks/Alias	RVP1; HRVP1; C7orf1; CPE-R2; CPETR2; Claudin-3
Species	Human
GeneID (Human)	1365
GeneID	1365
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 121-220 of human Claudin-3 mAb (NP_001297.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O15551
Protein Weight	23kDa
Shipping	Ice bag