

Anti-Claudin-4 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - O14493

Item Number: CM0005513



Introduction

CamelBio offers a high-quality rabbit monoclonal antibody against human Claudin-4, validated for WB, IP, and ELISA. Ideal for tight junction research and CPE receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Claudin-4 Rabbit mAb
Remarks/Alias	CLDN4; CPE-R; CPER; CPETR; CPETR1; WBSCR8; hCPE-R; claudin-4
Species	Human
GeneID (human)	1364
Immunogen	Synthetic peptide corresponding to amino acids 111-209 of human Claudin-4 (NP_001296.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	O14493
Protein Weight	22kDa
Shipping	Ice bag