

Anti-Claudin 5 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Ip, Elisa - O00501

Item Number: CM0011182

Introduction

Claudin 5 Rabbit mAb validated for WB, IF/ICC, IHC-P, IP, ELISA with cross-reactivity to human, mouse, rat. Ideal for tight junction research and cell junction studies.

[Learn More](#)

Parameter	Value
Product Name	Claudin 5 Rabbit mAb
Remarks/Alias	AWAL; BEC1; TMVCF; TMDVCF; CPETRL1
Species	Human
GeneID (Human)	7122
GeneID	7122
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 119-218 of human Claudin 5 (NP_001349995.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00501
Protein Weight	23kDa
Shipping	Ice bag