

Anti-Claudin 6 Polyclonal Antibody For Wb And Elisa - P56747

Item Number: CM0021016

Introduction

Rabbit polyclonal antibody targeting human Claudin 6 (CLDN6), validated for WB and ELISA. Cross-reacts with human, mouse, rat. Useful for tight junction and HCV entry studies in liver cells.

[Learn More](#)

Parameter	Value
Product Name	Claudin 6 Rabbit pAb
Remarks/Alias	CLDN6
Species	Human
Gene ID (Human)	9074
Gene ID	9074
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 121-220 of human Claudin 6 (NP_067018.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56747
Protein Weight	23kDa
Shipping	Ice bag