

Anti-Claudin 6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P56747

Item Number: CM0006448

Introduction

Rabbit monoclonal antibody targeting human Claudin 6. Suitable for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Essential for research on tight junctions and HCV entry.

[Learn More](#)

Specification	Detail
Product Name	Claudin 6 Rabbit mAb
Remarks/Alias	Skullin
Species	Human
GeneID (Human)	9074
GeneID	9074
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 119-219 of human Claudin 6 (NP_067018.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56747
Protein Weight	23kDa
Shipping	Ice bag