

## Abflo 647 Anti-Claudin18.2 Monoclonal Antibody For Flow Cytometry - P56856-2

Item Number: CM0009760

### Introduction

High-quality ABflo 647-conjugated rabbit monoclonal antibody targeting human Claudin18.2 (isoform A2). Validated for flow cytometry, recognizing a 27kDa tight junction protein expressed in gastric mucosa. Ideal for gastric barrier research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human Claudin18.2 mAb                                                                                         |
| Remarks/Alias    | SFTA5                                                                                                                     |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 51208                                                                                                                     |
| Gene ID          | 51208                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Claudin18.2 (NP_001002026.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | FC                                                                                                                        |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | P56856-2                                                                                                                  |
| Protein Weight   | 27kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |