


Anti-Clca1 Polyclonal Antibody For Wb, If-P, Elisa - A8K7I4

Item Number: CM0012875



Introduction

Rabbit polyclonal antibody targeting CLCA1 (Calcium-activated chloride channel regulator 1). Validated for WB, IF-P, ELISA; cross-reacts with human, mouse, rat. Ideal for mucus hypersecretion, inflammation, and ion transport studies.

[Learn More](#)

Parameter	Value
Product Name	CLCA1 Rabbit pAb
Remarks/Alias	CACC, GOB5, CACC1, CLCRG1, CaCC-1, hCLCA1, hCaCC-1, CLCA1
Species	Human
GeneID (Human)	1179
GeneID	1179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-150 of human CLCA1 (NP_001276.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A8K7I4
Protein Weight	100kDa
Shipping	Ice bag