

Anti-Clcc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96S66

Item Number: CM0021540

Introduction

Rabbit polyclonal antibody against human CLCC1, validated for WB, IHC-P, ELISA. Targets the ER chloride channel important for calcium homeostasis and neuronal survival. Cross-reacts with mouse. Suitable for ER stress and retinal research.

[Learn More](#)

Parameter	Value
Product Name	CLCC1 Rabbit pAb
Remarks/Alias	MCLC; RP32; CLCC1
Species	Human
GeneID (Human)	23155
GeneID	23155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-501 of human CLCC1 (NP_055942.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96S66
Protein Weight	62kDa
Shipping	Ice bag