


Anti-Clcn1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P35523

Item Number: CM0025612



Introduction

CLCN1 Rabbit pAb (P35523) is a polyclonal antibody targeting human chloride channel protein 1. Validated for WB, IHC-P, ELISA, with cross-reactivity to human, mouse, and rat. Supports skeletal muscle excitability and channelopathy research.

[Learn More](#)

Parameter	Value
Product Name	CLCN1 Rabbit pAb
Remarks/Alias	CLC1; CLCN1
Species	Human
GeneID (Human)	1180
GeneID	1180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 779-988 of human CLCN1 (NP_000074.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35523
Protein Weight	109kDa
Shipping	Ice bag