


Anti-Clcn2 Rabbit Polyclonal Antibody For Wb And Elisa - P51788

Item Number: CM0030369



Introduction

High-quality rabbit polyclonal antibody against human CLCN2, validated for WB and ELISA. Detects human and mouse chloride channel protein 2. Ideal for neuroscience, ion channel, and adrenal research.

[Learn More](#)

Parameter	Value
Product Name	CLCN2 Rabbit pAb
Remarks/Alias	CLC2; ECA2; ECA3; EGI3; EGMA; EJM6; EJM8; FHA2; FHII; CIC-2; EGI11; HALD2; LKPAT; cIC-2; CLCN2
Species	Human
Gene ID (Human)	1181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-770 of human CLCN2 (NP_001164560.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51788
Protein Weight	99kDa
Shipping	Ice bag