


## Anti-Clcn7 Rabbit Polyclonal Antibody For Wb And Elisa - P51798

Item Number: CM0027875



### Introduction

Anti-CLCN7 rabbit polyclonal antibody validated for Western blot and ELISA applications. Detects CLCN7 (CLC-7) protein in human, mouse, and rat samples. Suitable for lysosomal acidification, chloride transport, and osteopetrosis research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CLCN7 Rabbit pAb                                                                                                   |
| Remarks/Alias    | HOD; CLC7; CLC-7; OPTA2; OPTB4; PPP1R63; CLCN7                                                                     |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 1186                                                                                                               |
| GeneID           | 1186                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 626-805 of human CLCN7 (NP_001278.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P51798                                                                                                             |
| Protein Weight   | 89kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |