

Anti-Clon10 Rabbit Polyclonal Antibody For Wb, Elisa - P78369

Item Number: CM0027357

Introduction

Anti-CLDN10 rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to rat. Targets human Claudin-10, a tight junction protein that regulates paracellular ion permeability, suitable for renal and epithelial research.

[Learn More](#)

Parameter	Value
Product Name	CLDN10 Rabbit pAb
Remarks/Alias	OSPL; HELIX; OSP-L; CPETRL3; CLDN10
Species	Human
GeneID (Human)	9071
GeneID	9071
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLDN10 (NP_008915.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P78369
Protein Weight	24kDa
Shipping	Ice bag