

Anti-Clon11 Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - 075508

Item Number: CM0017791

Introduction

CLDN11 Rabbit Polyclonal Antibody designed for WB, IHC-P, ELISA applications. Recognizes human, mouse, and rat CLDN11, a tight junction component implicated in intercellular adhesion. Ideal for neurobiology and barrier function studies.

[Learn More](#)

Parameter	Value
Product Name	CLDN11 Rabbit pAb
Remarks/Alias	OSP; OTM; HLD22; CLDN11
Species	Human
Gene ID	5010 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLDN11 (NP_005593.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075508
Protein Weight	22kDa
Shipping	Ice bag