

Anti-ClDn17 Polyclonal Antibody For Wb, Ihc-P, Elisa - P56750

Item Number: CM0013221

Introduction

Rabbit polyclonal antibody targeting human Claudin-17 (CLDN17) for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for tight junction and anion transport research.

[Learn More](#)

Parameter	Specification
Product Name	CLDN17 Rabbit pAb
Remarks/Alias	CLDN17
Species	Human
Gene ID (Human)	26285
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human CLDN17 (NP_036263.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P56750
Protein Weight	25kDa
Shipping	Ice bag