

Anti-CLDN3 Rabbit Polyclonal Antibody For Wb, Elisa - 015551

Item Number: CM0016745

Introduction

Rabbit polyclonal antibody against human CLDN3, validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Useful for tight junction and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CLDN3 Rabbit pAb
Remarks/Alias	RVP1; HRVP1; C7orf1; CPE-R2; CPETR2; CLDN3
Species	Human
Gene ID	1365
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 141-220 of human CLDN3 (NP_001297.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	015551
Protein Weight	23kDa
Shipping	Ice bag