

Anti-Cldn5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00501

Item Number: CM0018336

Introduction

High-quality rabbit polyclonal antibody against CLDN5, suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Essential for tight junction and barrier function studies.

[Learn More](#)

Parameter	Value
Product Name	CLDN5 Rabbit pAb
Remarks/Alias	AWAL; BEC1; TMVCF; TMDVCF; CPETRL1; CLDN5
Species	Human
GeneID (Human)	7122
GeneID	7122
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 119-218 of human CLDN5 (NP_001349995.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00501
Protein Weight	23kDa
Shipping	Ice bag