

Anti-Cltn7 Polyclonal Antibody For Wb, If-P, Elisa - 095471

Item Number: CM0021855

Introduction

Anti-CLDN7 rabbit polyclonal antibody validated for WB, IF-P, and ELISA; detects human, mouse, and rat Claudin-7 at ~22 kDa. Ideal for tight junction studies and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	CLDN7 Rabbit pAb
Remarks/Alias	CLDN-7, CEPTRL2, CPETRL2, Hs.84359, claudin-1, CLDN7
Species	Human
Gene ID (Human)	1366
Gene ID	1366
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-211 of human CLDN7 (NP_001298.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95471
Protein Weight	22kDa
Shipping	Ice bag