


Anti-Cleaved Caspase-1 P20 Rabbit Monoclonal Antibody For Wb, Elisa - P29466

Item Number: CM0007082



Introduction

Rabbit monoclonal antibody targeting the cleaved p20 subunit of human caspase-1, validated for WB and ELISA. Suitable for studies of inflammasome activation, pyroptosis, and inflammatory cytokine processing.

[Learn More](#)

Parameter	Specification
Product Name	Cleaved Caspase-1 p20 Rabbit mAb
Remarks/Alias	ICE; P45; IL1BC; Caspase-1 p20
Species	Human
GeneID (Human)	834
GeneID	834
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Caspase-1 p20 (NP_001244047.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P29466
Protein Weight	10KDa/29KDa/35KDa/42KDa/45KDa
Shipping	Ice bag