


Anti-Cleaved-Caspase 3 P12 Rabbit Monoclonal Antibody For Wb, Elisa - P42574

Item Number: CM0011993



Introduction

High-specificity rabbit monoclonal antibody targeting the p12 cleavage fragment of human caspase-3 (P42574). Validated for WB and ELISA, supporting apoptosis mechanism studies and biomarker development.

[Learn More](#)

Parameter	Value
Product Name	Cleaved-Caspase 3 p12 Rabbit mAb
Remarks/Alias	CPP32; SCA-1; CPP32B
Species	Human
Gene ID (Human)	836
Gene ID	836
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 171-270 of human Cleaved-Caspase 3 p12 (NP_004337.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42574
Protein Weight	32 kDa
Shipping	Ice bag