


Anti-Cleaved Gasdermin D (N-Terminal) Rabbit Monoclonal Antibody For Pyroptosis Research - P57764

Item Number: CM0006306



Introduction

Highly specific rabbit monoclonal antibody targeting the cleaved N-terminal fragment of human Gasdermin D (P57764). Suitable for WB, IF/ICC, and ELISA. Ideal for studying pyroptosis, inflammasome activation, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cleaved Gasdermin D (N Terminal) Rabbit mAb
Remarks/Alias	DF5L; DFNA5L; FKSG10; GSDMDC1
Species	Human
Gene ID (Human)	79792
Gene ID	79792
Immunogen	Synthetic peptide corresponding to residues surrounding Asp275 of human Gasdermin D (NP_079012.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P57764
Protein Weight	53kDa
Shipping	Ice bag