


Anti-Cleaved Parp (Asp214) Monoclonal Antibody For Wb, If/Icc, Elisa - P09874

Item Number: CM0005559



Introduction

Rabbit monoclonal anti-Cleaved PARP (Asp214) antibody validated for WB, IF/ICC, ELISA. Detects human PARP1 apoptosis-specific cleavage fragment. Ideal for DNA damage and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Cleaved PARP (Asp214) Rabbit mAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP
Species	Human
Gene ID	142
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 201-300 of human PARP (Asp214) (NP_031441.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P09874
Protein Weight	113kDa
Shipping	Ice bag