

Anti-Clec11A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y240

Item Number: CM0025860

Introduction

Rabbit polyclonal antibody against human CLEC11A, validated for WB and ELISA, reacts with human and mouse. Ideal for studying osteogenesis, bone maintenance, and skeletal tissue research.

[Learn More](#)

Parameter	Value
Product Name	CLEC11A Rabbit pAb
Remarks/Alias	P47; SCGF; LSLCL; CLECSF3; CLEC11A
Species	Human
GeneID (human)	6320
GeneID	6320
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-323 of human CLEC11A (NP_002966.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y240
Protein Weight	36kDa
Shipping	Ice bag